

Surgical scrub sink MP'DUO® Electronic

Aseptic surgical scrub sink molded in one piece without any join or sealing.

With no overflow hole in order to avoid bacterial development, MP'DUO® has a large wall and sides protections against splashing water.

Available with a full range of accessories for an ideal configuration matching your specific need. For high risk services and critical medical units, this surgical sink is highly recommended with our filtering option for pathogen free water

Our team advises every customer according to the country's water quality and environment.

Material: Made of reinforced molded polyester resin fiberglass. The outer coating is made of high grade iso NPG sanitary gelcoat.

Targeted services: Operating rooms, emergency boxes, intensive care unit, dental surgery offices, cosmetic and plastic surgery offices, nephrology and dialysis units, critical care, neonatal unit.

Water command

Electronic infrared sensor




With anti-bacterial automatic flush (6h/12h/Off after last use)



Dimensions

L : 1500 mm
H : 745 mm
D : 570 mm
Weight : 19 Kg



Standard equipment

- 2 Infrared sensors with adjustable detection range
 - + 220V AC device and Battery device. 2 in 1
 - + Easy setting with an optional remote control
 - + Control box easy to unclip for maintenance
- 2 Gooseneck taps, chromed finish
- 1 Siphon of Ø40mm
- 2 Regulation valves
- HDPE flexible connecting hoses
- Stainless steel fixing frame

Optional accessories

- Straight faucet taps, premium heavy chromed finish
- Thermostatic mixer with anti-burn security
- Prefiltering kit 1µ with pressure protection system
- Ultraviolet kit for water disinfection - 40W
- Medical disinfectable glass siphon
- Terminal filters 0.2µ all germs 62 days
- Aluminum rail for accessories + 2 clipping frames
- Elbow command Airless Soap dispenser
- Electronic command Airless Soap dispenser
- Manual paper dispenser
- Nail brush dispenser
- Basket for accessories

Surgical scrub sink MP'DUO® Knee control

Aseptic surgical scrub sink molded in one piece without any join or sealing.

With no overflow hole in order to avoid bacterial development, MP'DUO® has a large wall and sides protections against splashing water.

Available with a full range of accessories for an ideal configuration matching your specific need. For high risk services and critical medical units, this surgical sink is highly recommended with our filtering option for pathogen free water

Our team advises every customer according to the country's water quality and environment.


Material: Made of reinforced molded polyester resin fiberglass. The outer coating is made of high grade iso NPG sanitary gelcoat.

Targeted services: Operating rooms, emergency boxes, intensive care unit, dental surgery offices, cosmetic and plastic surgery offices, nephrology and dialysis units, critical care, neonatal unit.



Dimensions

L : 1500 mm
H : 745 mm
D : 570 mm
Weight : 19 Kg



Water command

Mechanical knee control
Temporized 10s



Standard equipment

- 2 Knee control sticks with tempomatic valves
- 2 Gooseneck taps, chromed finish
- 1 Siphon of Ø40mm
- 2 Regulation valves
- HDPE flexible connecting hoses
- Stainless steel fixing frame

Optional accessories

- Straight faucet taps, premium heavy chromed finish
- Thermostatic mixer with anti-burn security
- Prefiltering kit 1µ with pressure protection system
- Ultraviolet kit for water disinfection - 40W
- Medical disinfectable glass siphon
- Terminal filters 0.2µ all germs 62 days
- Aluminum rail for accessories + 2 clipping frames
- Elbow command Airless Soap dispenser
- Electronic command Airless Soap dispenser
- Manual paper dispenser
- Nail brush dispenser
- Basket for accessories

MEDICAL PROCESS S.A.S.

contact@medicalprocess.fr | www.medicalprocess.fr

Tel. 0033(0)9.83.98.43.56 - Fax 0033(0)9.81.72.93.56

311 av des carrières, Z.I la Plaine du Caire 3

13830 Roquefort la Bédoule - FRANCE



Surgical scrub sink MP'DUO® Elbow control

Aseptic surgical scrub sink molded in one piece without any join or sealing.

With no overflow hole in order to avoid bacterial development, MP'DUO® has a large wall and sides protections against splashing water.

Available with a full range of accessories for an ideal configuration matching your specific need. For high risk services and critical medical units, this surgical sink is highly recommended with our filtering option for pathogen free water

Our team advises every customer according to the country's water quality and environment.


Material: Made of reinforced molded polyester resin fiberglass. The outer coating is made of high grade iso NPG sanitary gelcoat.

Targeted services: Operating rooms, emergency boxes, intensive care unit, dental surgery offices, cosmetic and plastic surgery offices, nephrology and dialysis units, critical care, neonatal unit.



Dimensions

L : 1500 mm
H : 745 mm
D : 570 mm
Weight : 19 Kg



Water command

Mechanical elbow control
Premium heavy chromed finish



Allows mixing manually Hot and Cold water

Standard equipment

- 2 Elbow control taps with manual mixing sticks
- 1 Siphon of Ø40mm
- 2 Regulation valve
- HDPE flexible connecting hoses
- Stainless steel fixing frame

Optional accessories

- Prefiltering kit 1µ with pressure protection system
- Ultraviolet kit for water disinfection - 40W
- Medical disinfectable glass siphon
- Terminal filters 0.2µ all germs 62 days
- Aluminum rail for accessories + 2 clipping frames
- Elbow command Airless Soap dispenser
- Electronic command Airless Soap dispenser
- Manual paper dispenser
- Nail brush dispenser
- Basket for accessories

MEDICAL PROCESS S.A.S.

contact@medicalprocess.fr | www.medicalprocess.fr

Tel. 0033(0)9.83.98.43.56 - Fax 0033(0)9.81.72.93.56

311 av des carrières, Z.I la Plaine du Caire 3

13830 Roquefort la Bédoule - FRANCE

